

DIASWISS

Products for
cosmetic podiatry and manicure

Produits pour
cosmétiques podologie et manicure

Produkte für
die kosmetische Fuss- & Nagelpflege

Productos para el cuidado
cosmético de los pies y las uñas

Prodotti per
la cura cosmetica di piedi e unghie



Tel: +41 22 301 56 00 · Fax: +41 22 301 56 50

mail@diaswiss.ch · www.diaswiss.ch

Electroplated diamonds Instrumentos de diamante galvánicos

Instruments diamantés galvaniques Strumenti diamantati galvanizzati

Galvanische Diamanten

Shank **HP**
Speed 5' - 30.000 rpm



F	shape	D79	D79	D79	D104	D251	D801	D840	D848	D848	D860	D881
ISO	ISO-code	261	200	200	104	275	001	111	172	172	245	141
Ø 1/10 mm	head diameter	050	063	063	070	080	023	040	018	040	014	016
L /mm	head length	12	14	14	3	15	—	9	8	9	5	7
D	diamond coating	M	SG	M	M	M	M	M	M	M	M	M
P	price category	14	14	14	14	14	12	14	11	14	11	11

Order-No

D79 050
D79 SG 063
D79 063
D104 070
D251 080
D801 023
D840 040
D848 018
D848 040
D860 014
D881 016

Turbo diamonds with spiral Fresas diamantadas con espiral

Fraises diamantées avec spirale Frese diamantate a spirale

Turbodiamanten mit Spirale

Shank **HP**
Speed 3' - 15.000 rpm



F	shape	D144	D201	D266	D881	D881	D881
ISO	ISO-code	144	201	266	143	144	145
Ø 1/10 mm	head diameter	110	018	062	080	105	120
L /mm	head length	14	10	20	13	15	18
P	price category	101	12	101	101	101	101

Order-No

D144 110
D201 018
D266 062
D881 080
D881 105
D881 120

Carbide burs Fresas de carburo de tungsteno

Fraises en carbure de tungstène Frese in carburo di tungsteno

Hartmetallfräser

Shank **HP**
Speed 5' - 30.000 rpm



F	shape	D72	D251	D251	D251	D251	D253
ISO	ISO-code	144	275	275	292	275	277
Ø 1/10 mm	head diameter	060	060	060	060	060	014
L /mm	head length	13	10,5	10,5	14	10,5	3,5
P	price category	53	52	52	50	54	41

Order-No

D72 F 060
D251 C 060
D251 F 060
D251 G 060
D251 ST 060
D253 F 014

Packaging unit: 1
Unité d'emballage: 1
Verpackungseinheit: 1
Unidad de envasado: 1
Unità d'imballaggio: 1